

FNBP3 Antibody

Catalog No: #47954

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FNBP3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of FNBP3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human FNBP3.
Target Name	PRPF40A
Other Names	FBP11; FLAF1; FNBP3; HIP10; HYPA; Pre-mRNA-processing factor 40 homolog A; Fas ligand-associated factor 1; Formin-binding protein 11; Formin-binding protein 3; Huntingtin yeast partner A; Huntingtin-interacting protein 10; HIP-10; Huntingtin-interacting protein A; Renal carcinoma antigen NY-REN-6
Accession No.	Swiss-Prot#:O75400NCBI Gene ID:55660
Uniprot	O75400
GeneID	55660;
Calculated MW	108KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of FNBP3 expression in HeLa (A), mouse brain (B), rat brain (C) whole cell lysates.

Note: This product is for in vitro research use only