

NIF1 Antibody

Catalog No: #47955

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NIF1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of NIF1 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human NIF1.
Target Name	ZNF335
Other Names	Zinc finger protein 335; NRC-interacting factor 1; NIF-1
Accession No.	Swiss-Prot#:Q9H4Z2NCBI Gene ID:63925
Uniprot	Q9H4Z2
GeneID	63925;
Calculated MW	144KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of NIF1 expression in HEK293T (A), RAW264.7 (B), H9C2 (C) whole cell lysates

Note: This product is for in vitro research use only