

TBX2/3 Antibody

Catalog No: #47959

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TBX2/3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of TBX2/3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human TBX2/3.
Target Name	TBX2; TBX3
Other Names	TBX3; T-box transcription factor TBX3; T-box protein 3; TBX2; T-box transcription factor TBX2; T-box protein 2
Accession No.	Swiss-Prot#:O15119; Q13207NCBI Gene ID:6926
Uniprot	O15119
GeneID	6926;
Calculated MW	79KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of TBX3 expression in CT26 (A), MEF (B), C6 (C), HEK293T (D), HepG2 (E) whole cell lysates

Note: This product is for in vitro research use only