

## Claudin 17 Antibody

Catalog No: #47960

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Claudin 17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of Claudin 17 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Claudin 17.
Target Name	CLDN17
Other Names	Claudin-17
Accession No.	Swiss-Prot#:P56750NCBI Gene ID:26285
Uniprot	P56750
GeneID	26285;
Calculated MW	25KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/2000)

## Images



Western blot analysis of Claudin 17 expression in KEK293T (A), PC12 (B), MEF (C) whole cell lysates

Note: This product is for in vitro research use only