

ZNF307 Antibody

Catalog No: #47961

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF307 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IHC
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of ZNF307 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ZNF307.
Target Name	ZKSCAN4
Other Names	ZNF307; ZNF427; Zinc finger protein with KRAB and SCAN domains 4; P373c6.1; Zinc finger protein 307; Zinc finger protein 427
Accession No.	Swiss-Prot#:Q969J2NCBI Gene ID:387032
Uniprot	Q969J2
GeneID	387032;
Calculated MW	65KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C


Application Details

WB (1/500 - 1/2000), IHC (1/50 - 1/200)

Images



Immunohistochemical analysis of ZNF307 staining in human brain formalin fixed paraffin embedded tissue section. The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



Western blot analysis of ZNF307 expression in MEF (A), PC12 (B), A549 (C) whole cell lysates

Note: This product is for in vitro research use only