

## Ribosomal Protein L26L Polyclonal Antibody

Catalog No: #41403

Package Size: #41403-1 50ul #41403-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Ribosomal Protein L26L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	Ribosomal Protein L26L Polyclonal Antibody detects endogenous levels of Ribosomal Protein L26L protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L26L.
Target Name	Ribosomal Protein L26L
Other Names	RPL26L1; RPL26P1; 60S ribosomal protein L26-like 1
Accession No.	Swiss-Prot: Q9UNX3NCBI Gene ID: 51121
Uniprot	Q9UNX3
GeneID	51121;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

## Images



Western Blot analysis of MCF7 cells using Ribosomal Protein L26L Polyclonal Antibody

---

Note: This product is for in vitro research use only