

TWIK3 Antibody

Catalog No: #47968

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | TWIK3 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Applications | WB |
| Species Reactivity | Hu,Ms,Rt |
| Specificity | Recognizes endogenous levels of TWIK3 protein. |
| Immunogen Description | KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human TWIK3. |
| Target Name | KCNK7 |
| Other Names | Potassium channel subfamily K member 7 |
| Accession No. | Swiss-Prot#:Q9Y2U2NCBI Gene ID:10089 |
| Uniprot | Q9Y2U2 |
| GeneID | 10089; |
| Calculated MW | 36KD |
| Concentration | 1 mg/ml |
| Formulation | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage | Store at -20°C |

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of TWIK3 expression in HCT116 (A), A549 (B), H9C2 (C), AML12 (D) whole cell lysates

Note: This product is for in vitro research use only