

VPS26B Antibody

Catalog No: #47971

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VPS26B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of VPS26B protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human VPS26B.
Target Name	VPS26B
Other Names	Vacuolar protein sorting-associated protein 26B; Vesicle protein sorting 26B
Accession No.	Swiss-Prot#:Q4G0F5NCBI Gene ID:112936
Uniprot	Q4G0F5
GeneID	112936;
Calculated MW	39KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of VPS26B expression in MCF7 (A), HeLa (B), PC12 (C), AML12 (D) whole cell lysates.

Note: This product is for in vitro research use only