

hnRNP CL1 Antibody

Catalog No: #47975

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	hnRNP CL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms
Specificity	Recognizes endogenous levels of hnRNP CL1 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human hnRNP CL1.
Target Name	HNRNPCL1
Other Names	HNRNPCL1; Heterogeneous nuclear ribonucleoprotein C-like 1; hnRNP C-like-1; hnRNP core protein C-like 1
Accession No.	Swiss-Prot#:O60812NCBI Gene ID:343069
Uniprot	O60812
GeneID	343069;
Calculated MW	32KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of hnRNP CL1 expression in A549 (A), HEK293T (B), MEF (C) whole cell lysates

Note: This product is for in vitro research use only