

Kv5.1 Antibody

Catalog No: #47977

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Kv5.1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of Kv5.1 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Kv5.1.
Target Name	KCNF1
Other Names	Potassium voltage-gated channel subfamily F member 1; Voltage-gated potassium channel subunit Kv5.1; kH1
Accession No.	Swiss-Prot#:Q9H3M0NCBI Gene ID:3754
Uniprot	Q9H3M0
GeneID	3754;
Calculated MW	55KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of Kv5.1 expression in COS7 (A), Panc1 (B), mouse heart (C), rat heart (D) whole cell lysates.

Note: This product is for in vitro research use only