

# FAM35A/B Antibody

Catalog No: #47978



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FAM35A/B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of FAM35A/B protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human FAM35A/B.
Target Name	FAM35A; FAM35B
Other Names	FAM35A; Protein FAM35A; FAM35B; Protein FAM35B
Accession No.	Swiss-Prot#:Q86V20; A6NDK8NCBI Gene ID:54537
Uniprot	Q86V20
GeneID	54537;
Calculated MW	94KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of FAM35A/B expression in H1792 (A), H9C2 (B), AML12 (C) whole cell lysates.

Note: This product is for in vitro research use only