

RPAP2 Antibody

Catalog No: #47981

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RPAP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms
Specificity	Recognizes endogenous levels of RPAP2 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human RPAP2.
Target Name	RPAP2
Other Names	C1orf82; Putative RNA polymerase II subunit B1 CTD phosphatase RPAP2; RNA polymerase II-associated protein 2
Accession No.	Swiss-Prot#:Q8IXW5NCBI Gene ID:79871
Uniprot	Q8IXW5
GeneID	79871;
Calculated MW	70KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of RPAP2 expression in HEK293T (A), LO2 (B), A549 (C), AML12 (D), BV2 (E) whole cell lysates.

Note: This product is for in vitro research use only