

PDIK1L Antibody

Catalog No: #47982

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PDIK1L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB, IHC
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of PDIK1L protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PDIK1L.
Target Name	PDIK1L
Other Names	CLIK1L; Serine/threonine-protein kinase PDIK1L; PDLIM1-interacting kinase 1-like
Accession No.	Swiss-Prot#:Q8N165NCBI Gene ID:149420
Uniprot	Q8N165
GeneID	149420;
Calculated MW	39KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C


Application Details

WB (1/500 - 1/1000), IHC (1/50 - 1/100)

Images



Western blot analysis of PDIK1L expression in LO2 (A), HEK293T (B), COS7 (C) whole cell lysates



Immunohistochemical analysis of PDIK1L staining in human breast cancer formalin fixed paraffin embedded tissue section. The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Note: This product is for in vitro research use only