

Ribosomal Protein L28 Polyclonal Antibody

Catalog No: #41405

Package Size: #41405-1 50ul #41405-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ribosomal Protein L28 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Ribosomal Protein L28 Polyclonal Antibody detects endogenous levels of Ribosomal Protein L28 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ribosomal Protein L28.
Target Name	Ribosomal Protein L28
Other Names	RPL28; 60S ribosomal protein L28
Accession No.	Swiss-Prot: P46779NCBI Gene ID: 6158
Uniprot	P46779
GeneID	6158;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Ribosomal Protein L28 Polyclonal Antibody

Note: This product is for in vitro research use only