

Serpin B12 Antibody

Catalog No: #47983

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Serpin B12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms
Specificity	Recognizes endogenous levels of Serpin B12 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Serpin B12.
Target Name	SERPINB12
Other Names	Serpin B12
Accession No.	Swiss-Prot#:Q96P63NCBI Gene ID:89777
Uniprot	Q96P63
GeneID	89777;
Calculated MW	46KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of Serpin B12 expression in A549 (A), LO2 (B), AML12 (C) whole cell lysates.

Note: This product is for in vitro research use only