

CD172g Antibody

Catalog No: #47985

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD172g Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of CD172g protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human CD172g.
Target Name	SIRPG
Other Names	SIRPB2; Signal-regulatory protein gamma; SIRP-gamma; CD172 antigen-like family member B; Signal-regulatory protein beta-2; SIRP-b2; SIRP-beta-2; CD172g
Accession No.	Swiss-Prot#:Q9P1W8NCBI Gene ID:55423
Uniprot	Q9P1W8
GeneID	55423;
Calculated MW	42KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of CD172g expression in AML12 (A), PC12 (B), LO2 (C), HEK293T (D) whole cell lysates.

Note: This product is for in vitro research use only