

## GPR103 Antibody

Catalog No: #47988

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GPR103 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of GPR103 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human GPR103.
Target Name	QRFP
Other Names	GPR103; Pyroglutamylated RFamide peptide receptor; AQ27; G-protein coupled receptor 103; Orexigenic neuropeptide QRFP receptor; SP9155
Accession No.	Swiss-Prot#:Q96P65NCBI Gene ID:84109
Uniprot	Q96P65
GeneID	84109;
Calculated MW	49KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of GPR103 expression in U87MG (A), HEK293T (B) whole cell lysates.

Note: This product is for in vitro research use only