

Ribosomal Protein L7 Polyclonal Antibody

Catalog No: #41408

Package Size: #41408-1 50ul #41408-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribosomal Protein L7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Ribosomal Protein L7 Polyclonal Antibody detects endogenous levels of Ribosomal Protein L7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L7.
Target Name	Ribosomal Protein L7
Other Names	RPL7; 60S ribosomal protein L7
Accession No.	Swiss-Prot: P18124NCBI Gene ID: 6129
Uniprot	P18124
GeneID	6129;
SDS-PAGE MW	32kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Ribosomal Protein L7 Polyclonal Antibody

Note: This product is for in vitro research use only