

Ribosomal Protein S8 Polyclonal Antibody

Catalog No: #41411

Package Size: #41411-1 50ul #41411-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Ribosomal Protein S8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	Ribosomal Protein S8 Polyclonal Antibody detects endogenous levels of Ribosomal Protein S8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ribosomal Protein S8.
Target Name	Ribosomal Protein S8
Other Names	RPS8; OK/SW-cl.83; 40S ribosomal protein S8
Accession No.	Swiss-Prot: P62241NCBI Gene ID: 6202
Uniprot	P62241
GeneID	6202;
SDS-PAGE MW	30kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of COLO205 cells using Ribosomal Protein S8 Polyclonal Antibody

Note: This product is for in vitro research use only