

## V5-Tag Antibody

Catalog No: #48047

Package Size: #48047-1 50ul #48047-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	V5-Tag Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Immunogen Description	peptide
Other Names	GKPIPPELLGLDST epitope tag antibody GKPIPPELLGLDST tag antibody Protein Rev antibody Regulator of expression of viral proteins antibody rev antibody V5 epitope tag antibody
Accession No.	Swiss-Prot#:
Calculated MW	50 kDa
Formulation	1*TBS (pH7.4), 0.5%BSA, 25%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

## Application Details

WB: 1:10,000:100,000

## Images



Western blot analysis on V5-tagged recombinant protein (50kDa) using anti-V5 polyclonal antibody.

## Background

The V5 epitope tag is derived from a small epitope (PK) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The anti-V5 antibody detects recombinant proteins containing the 14 amino acid GKPIPPELLGLDST sequence.

## References

---

Note: This product is for in vitro research use only