

BSA Antibody

Catalog No: #48061

Package Size: #48061-1 50ul #48061-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BSA Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2-B8
Purification	ProA affinity purified
Applications	WB
Species Reactivity	Bovine
Immunogen Description	protein
Accession No.	Swiss-Prot#:P02769
Uniprot	P02769
GeneID	280717;
Calculated MW	66kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:500-1:1,000

Images



Western blot analysis on BSA using anti-BSA Mouse mAb (Cat. # M0809-3).

Background

Albumin is a soluble, monomeric protein has molecular weight of c. 66 kDa. It is primarily as a carrier protein for steroids, fatty acids, and thyroid hormones and plays a role in stabilizing extracellular fluid volume. Serum albumin protein has numerous biochemical applications including ELISAs, immunoblots, and immunohistochemistry. It is also used as a nutrient in cell and microbial culture.

References

Note: This product is for in vitro research use only