

Transmembrane protein 200A Antibody

Catalog No: #48081



Package Size: #48081-1 50ul #48081-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transmembrane protein 200A Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4-B3
Purification	Protein affinity purified
Applications	WB, ICC, IHC
Species Reactivity	Hu, Ms
Immunogen Description	Recombinant protein within human Transmembrane protein 200A 50-400 aa.
Other Names	KIAA1913 antibody T200A_HUMAN antibody Tmem200a antibody Transmembrane protein 200A antibody TTMC antibody Two transmembrane C antibody
Accession No.	Swiss-Prot#:Q86VY9
Uniprot	Q86VY9
GeneID	114801;
Calculated MW	47 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:500 IHC: 1:50 ICC: 1:100

Background

Transmembrane protein 200A(TM200A) is multi-pass membrane protein, which belongs to the TM200 family. It is expressed in cerebellum.

References

Note: This product is for in vitro research use only