

C1orf175 Antibody

Catalog No: #48086

Package Size: #48086-1 50ul #48086-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C1orf175 Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	A5-D4-A6
Purification	ProA affinity purified
Applications	WB, IHC
Species Reactivity	Hu
Immunogen Description	Recombinant protein
Accession No.	Swiss-Prot#:Q68CQ1
Uniprot	Q68CQ1
GeneID	374977;
Calculated MW	145 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:500-1:1000IHC: 1:50-1:200

Background

C1orf175 belongs to MROH family. MROH7 (Maestro Heat Like Repeat Family Member 7) is a Protein Coding gene. GO annotations related to this gene include binding. An important paralog of this gene is MROH7-TTC4.

References

Note: This product is for in vitro research use only