

RRS1 Polyclonal Antibody

Catalog No: #41417

Package Size: #41417-1 50ul #41417-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RRS1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	RRS1 Polyclonal Antibody detects endogenous levels of RRS1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human RRS1.
Target Name	RRS1
Other Names	RRS1; KIAA0112; RRR; Ribosome biogenesis regulatory protein homolog
Accession No.	Swiss-Prot: Q15050NCBI Gene ID: 23212
Uniprot	Q15050
GeneID	23212;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using RRS1 Polyclonal Antibody

Note: This product is for in vitro research use only