

Rsk-4 Polyclonal Antibody

Catalog No: #41418

Package Size: #41418-1 50ul #41418-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rsk-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Rsk-4 Polyclonal Antibody detects endogenous levels of Rsk-4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Rsk-4.
Target Name	Rsk-4
Other Names	RPS6KA6; RSK4; Ribosomal protein S6 kinase alpha-6; S6K-alpha-6; 90 kDa ribosomal protein S6 kinase 6; p90-RSK 6; p90RSK6; Ribosomal S6 kinase 4; RSK-4; pp90RSK4
Accession No.	Swiss-Prot: Q9UK32NCBI Gene ID: 27330
Uniprot	Q9UK32
GeneID	27330;
SDS-PAGE MW	84kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of 453 293T AD293 HELA HepG2 cells using Rsk-4 Polyclonal Antibody



Note: This product is for in vitro research use only