

## DPY30 Antibody

Catalog No: #48106



Package Size: #48106-1 50ul #48106-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DPY30 Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	ProA affinity purified
Applications	WB, ICC
Species Reactivity	Hu, Ms
Immunogen Description	peptide
Other Names	Dpy-30-like protein antibody Dpy-30L antibody DPY30 antibody DPY30_HUMAN antibody HDPY-30 antibody Protein dpy-30 homolog antibody Saf19 antibody
Accession No.	Swiss-Prot#:Q9C005
Uniprot	Q9C005
GeneID	84661;
Calculated MW	37kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

## Application Details

WB: 1:500

## Background

DPPA2 (developmental pluripotency associated 2) is an embryo-cancer antigen that is coexpressed with cancer-testis antigens in non-small cell lung cancer. Dppa2 is closely linked SAP motif genes restricted to pluripotent cells and the germ line. It may play a role in maintaining cell pluripotentiality and in embryo preimplantation development.

## References

Note: This product is for in vitro research use only