

# 6X His-Tag HRP Antibody

Catalog No: #48140

Package Size: #48140-1 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	6X His-Tag HRP Antibody
Purification	Immunogen affinity purified
Applications	WB
Immunogen Description	peptide
Conjugates	HRP-conjugated
Accession No.	Swiss-Prot#:
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

## Application Details

WB: 1:2,000-1:5,000

## Images




Fig1 Western blot analysis of extracts from cells expressing C-terminal His-tagged protein 6X His-tag antibody: Lane1: 1:2,000 diluted antibody Lane2: 1:5,000 diluted antibody

## Background

HRP conjugated His-Tag Antibody detects recombinant proteins containing the 6xHis epitope tag. The antibody recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.

Note: This product is for in vitro research use only