

Septin 8 Polyclonal Antibody

Catalog No: #41427

Package Size: #41427-1 50ul #41427-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Septin 8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Septin 8 Polyclonal Antibody detects endogenous levels of Septin 8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Septin 8.
Target Name	Septin 8
Other Names	SEPT8; KIAA0202; Septin-8
Accession No.	Swiss-Prot: Q92599NCBI Gene ID: 23176
Uniprot	Q92599
GeneID	23176;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Septin 8 Polyclonal Antibody

Note: This product is for in vitro research use only