

SLU7 Polyclonal Antibody

Catalog No: #41440

Package Size: #41440-1 50ul #41440-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLU7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SLU7 Polyclonal Antibody detects endogenous levels of SLU7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human SLU7.
Target Name	SLU7
Other Names	SLU7; Pre-mRNA-splicing factor SLU7; hSlu7
Accession No.	Swiss-Prot: O95391NCBI Gene ID: 10569
Uniprot	O95391
GeneID	10569;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of SH-SY5Y cells using SLU7 Polyclonal Antibody

Note: This product is for in vitro research use only