

SMYD2 Polyclonal Antibody

Catalog No: #41448

Package Size: #41448-1 50ul #41448-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SMYD2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	SMYD2 Polyclonal Antibody detects endogenous levels of SMYD2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human SMYD2.
Target Name	SMYD2
Other Names	SMYD2; KMT3C; N-lysine methyltransferase SMYD2; HSKM-B; Histone methyltransferase SMYD2; Lysine N-methyltransferase 3C; SET and MYND domain-containing protein 2
Accession No.	Swiss-Prot: Q9NRG4NCBI Gene ID: 56950
Uniprot	Q9NRG4
GeneID	56950;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using SMYD2 Polyclonal Antibody

Note: This product is for in vitro research use only