

Spindlin-1 Polyclonal Antibody

Catalog No: #41455

Package Size: #41455-1 50ul #41455-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Spindlin-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu Ms Rt
Specificity	Spindlin-1 Polyclonal Antibody detects endogenous levels of Spindlin-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Spindlin-1.
Target Name	Spindlin-1
Other Names	SPIN1; OCR; SPIN; Spindlin-1; Ovarian cancer-related protein
Accession No.	Swiss-Prot: Q9Y657NCBI Gene ID: 10927
Uniprot	Q9Y657
GeneID	10927;
SDS-PAGE MW	30kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of HeLa cells using Spindlin-1 Polyclonal Antibody



Western Blot analysis of 293T cells using Spindlin-1 Polyclonal Antibody

Note: This product is for in vitro research use only