

SR-1F Polyclonal Antibody

Catalog No: #41457

Package Size: #41457-1 50ul #41457-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	SR-1F Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt Mk
Specificity	SR-1F Polyclonal Antibody detects endogenous levels of SR-1F protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human SR-1F.
Target Name	SR-1F
Other Names	HTR1F; HTR1EL; 5-hydroxytryptamine receptor 1F; 5-HT-1F; 5-HT1F; Serotonin receptor 1F
Accession No.	Swiss-Prot: P30939NCBI Gene ID: 3355
Uniprot	P30939
GeneID	3355;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of COS7 cells using SR-1F Polyclonal Antibody

Note: This product is for in vitro research use only