

SRF Polyclonal Antibody

Catalog No: #41459

Package Size: #41459-1 50ul #41459-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SRF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	SRF Polyclonal Antibody detects endogenous levels of SRF protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human SRF around the non-phosphorylation site of S77.
Target Name	SRF
Other Names	SRF; Serum response factor; SRF
Accession No.	Swiss-Prot: P11831NCBI Gene ID: 6722
Uniprot	P11831
GeneID	6722;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MOUSE-BRAIN COLO cells using SRF Polyclonal Antibody

Note: This product is for in vitro research use only