

## GFP11-tag Antibody

Catalog No: #48545

Package Size: #48545-1 50ul #48545-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GFP11-tag Antibody
Purification	Peptide affinity purified.
Applications	WB
Immunogen Description	Synthetic peptide RDHMLVLEHYVNAAGITC linked to KLH.
Concentration	1mg/ml
Formulation	1*PBS (pH7.4), 0.2% BSA, 50% Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Application Details

WB: 1:1,000

## Images



Western blot analysis of GFP11-tag on fusion protein using anti-GFP11-tag antibody at 1/1,000 dilution.

## Background

Tagging proteins with a functional sequence is a powerful way to access function, and inserting tags at endogenous genomic loci allows the preservation of a near-native cellular background. GFP11 is the 11th beta-strand of the superfolder GFP beta-barrel structure. When expressed in the same cell, GFP11 and its complementary GFP fragment (GFP1-10) enable functional GFP tagging upon complementation. The key advantage of the GFP11 is its small size (16 aa): this allows commercial ssDNA oligomers to be used as HDR donors, circumventing any requirement for molecular cloning.

Note: This product is for in vitro research use only