

TADA2L Polyclonal Antibody

Catalog No: #41476

Package Size: #41476-1 50ul #41476-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TADA2L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	TADA2L Polyclonal Antibody detects endogenous levels of TADA2L protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TADA2L.
Target Name	TADA2L
Other Names	TADA2A; TADA2L; KL04P; Transcriptional adapter 2-alpha; Transcriptional adapter 2-like; ADA2-like protein
Accession No.	Swiss-Prot: O75478NCBI Gene ID: 6871
Uniprot	O75478
GeneID	6871;
SDS-PAGE MW	47kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using TADA2L Polyclonal Antibody

Note: This product is for in vitro research use only