

Tetraspanin-8 Polyclonal Antibody

Catalog No: #41487

Package Size: #41487-1 50ul #41487-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tetraspanin-8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Tetraspanin-8 Polyclonal Antibody detects endogenous levels of Tetraspanin-8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Tetraspanin-8.
Target Name	Tetraspanin-8
Other Names	TSPAN8; TM4SF3; Tetraspanin-8; Tspan-8; Transmembrane 4 superfamily member 3; Tumor-associated antigen CO-029
Accession No.	Swiss-Prot: P19075NCBI Gene ID: 7103
Uniprot	P19075
GeneID	7103;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of COLO205 cells using Tetraspanin-8 Polyclonal Antibody



Note: This product is for in vitro research use only