

TFIID Polyclonal Antibody

Catalog No: #41489

Package Size: #41489-1 50ul #41489-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TFIID Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	TFIID Polyclonal Antibody detects endogenous levels of TFIID protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TFIID.
Target Name	TFIID
Other Names	TBP; GTF2D1; TF2D; TFIID; TATA-box-binding protein; TATA sequence-binding protein; TATA-binding factor; TATA-box factor; Transcription initiation factor TFIID TBP subunit
Accession No.	Swiss-Prot: P20226NCBI Gene ID: 6908
Uniprot	P20226
GeneID	6908;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of 293 cells using TFIID Polyclonal Antibody



Note: This product is for in vitro research use only