

TGF α Polyclonal Antibody

Catalog No: #41490

Package Size: #41490-1 50ul #41490-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGF α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	TGF α Polyclonal Antibody detects endogenous levels of TGF α protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human TGF α .
Target Name	TGF α
Other Names	TGFA; Protransforming growth factor alpha
Accession No.	Swiss-Prot: P01135NCBI Gene ID: 7039
Uniprot	P01135
GeneID	7039;
SDS-PAGE MW	10kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of COS7 cells using TGF α Polyclonal Antibody

Note: This product is for in vitro research use only