

mTOR Rabbit mAb

Catalog No: #48834

Package Size: #48834-1 50ul #48834-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	mTOR Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	SU30-00
Purification	ProA affinity purified
Applications	WB, IHC, IP
Species Reactivity	Hu, Ms, Rt
Immunogen Description	recombinant protein
Other Names	dJ576K7.1 (FK506 binding protein 12 rapamycin associated protein 1) antibody FK506 binding protein 12 rapamycin associated protein 1 antibody FK506 binding protein 12 rapamycin associated protein 2 antibody FK506 binding protein 12 rapamycin complex associated protein 1 antibody FK506-binding protein 12-rapamycin complex-associated protein 1 antibody FKBP rapamycin associated protein antibody FKBP12 rapamycin complex associated protein antibody FKBP12-rapamycin complex-associated protein 1 antibody FKBP12-rapamycin complex-associated protein antibody FLJ44809 antibody FRAP antibody FRAP1 antibody FRAP2 antibody Mammalian target of rapamycin antibody Mechanistic target of rapamycin antibody mTOR antibody MTOR_HUMAN antibody OTTHUMP0000001983 antibody RAFT1 antibody Rapamycin and FKBP12 target 1 antibody Rapamycin associated protein FRAP2 antibody Rapamycin target protein 1 antibody Rapamycin target protein antibody RAPT1 antibody Serine/threonine-protein kinase mTOR antibody
Accession No.	Swiss-Prot#:P42345
Uniprot	P42345
GeneID	2475;
Calculated MW	289 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C


Application Details

WB: 1:1,000 IHC: 1:50-1:200


Images




Western blot analysis of mTOR on mouse testis lysates using anti-mTOR antibody at 1/1,000 dilution.




Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti-mTOR antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-mTOR antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue using anti-mTOR antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue using anti-mTOR antibody. Counter stained with hematoxylin.

Background

The phosphatidylinositol kinase (PIK) family members fall into two distinct subgroups. The first subgroup contains proteins such as the PI 3- and PI 4-kinases and the second group comprises the PIK-related kinases. The PIK-related kinases include Atm, DNA-PKCS and FRAP. These proteins have in common a region of homology at their carboxy-termini that is not present in the PI 3- and PI 4-kinases. The Atm gene is mutated in the autosomal recessive disorder ataxia telangiectasia (AT) that is characterized by cerebellar degeneration (ataxia) and the appearance of dilated blood vessels (telangiectases) in the conjunctivae of the eyes. AT cells are hypersensitive to ionizing radiation, impaired in mediating the inhibition of DNA synthesis and display delays in p53 induction. DNA-PK is a heterotrimeric DNA binding enzyme that is composed of a large subunit, DNA-PKCS, and two smaller subunits collectively known as Ku. The loss of DNA-PK leads to defects in DSB repair and V(D)J recombination. FRAP can autophosphorylate on serine and bind to Rapamycin/FKBP. FRAP is also an upstream regulator of S6 kinase and has been implicated in the regulation of p27 and p21 expression.

References

Note: This product is for in vitro research use only