

TGF β RII Polyclonal Antibody

Catalog No: #41493

Package Size: #41493-1 50ul #41493-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGF β RII Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	TGF β RII Polyclonal Antibody detects endogenous levels of TGF β RII protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TGF β RII.
Target Name	TGF β RII
Other Names	TGFBR2; TGF-beta receptor type-2; TGFR-2; TGF-beta type II receptor; Transforming growth factor-beta receptor type II; TGF-beta receptor type II; TbetaR-II
Accession No.	Swiss-Prot: P37173NCBI Gene ID: 7048
Uniprot	P37173
GeneID	7048;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using TGFβ RII Polyclonal Antibody

Note: This product is for in vitro research use only