

## TGFβ2 Polyclonal Antibody

Catalog No: #41495

Package Size: #41495-1 50ul #41495-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TGFβ2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	TGFβ2 Polyclonal Antibody detects endogenous levels of TGFβ2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human TGFβ2.
Target Name	TGFβ2
Other Names	TGFB2; Transforming growth factor beta-2; TGF-beta-2; BSC-1 cell growth inhibitor; Cetermin; Glioblastoma-derived T-cell suppressor factor; G-TSF; Polyergin
Accession No.	Swiss-Prot: P61812NCBI Gene ID: 7042
Uniprot	P61812
GeneID	7042;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

## Images

Western Blot analysis of mouse-brain cells using TGFβ2 Polyclonal Antibody



Note: This product is for in vitro research use only