

TGFβ3 Polyclonal Antibody

Catalog No: #41496

Package Size: #41496-1 50ul #41496-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TGFβ3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	TGFβ3 Polyclonal Antibody detects endogenous levels of TGFβ3 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human TGFβ3.
Target Name	TGFβ3
Other Names	TGFB3; Transforming growth factor beta-3; TGF-beta-3
Accession No.	Swiss-Prot: P10600NCBI Gene ID: 7043
Uniprot	P10600
GeneID	7043;
SDS-PAGE MW	13kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using TGFβ3 Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using TGF beta3 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only