

TIMP-3 Polyclonal Antibody

Catalog No: #41501

Package Size: #41501-1 50ul #41501-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TIMP-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	TIMP-3 Polyclonal Antibody detects endogenous levels of TIMP-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TIMP-3.
Target Name	TIMP-3
Other Names	TIMP3; Metalloproteinase inhibitor 3; Protein MIG-5; Tissue inhibitor of metalloproteinases 3; TIMP-3
Accession No.	Swiss-Prot: P35625NCBI Gene ID: 7078
Uniprot	P35625
GeneID	7078;
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using TIMP-3 Polyclonal Antibody

Note: This product is for in vitro research use only