

TNF- α Polyclonal Antibody

Catalog No: #41504

Package Size: #41504-1 50ul #41504-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF- α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	TNF- α Polyclonal Antibody detects endogenous levels of TNF- α protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TNF- α .
Target Name	TNF- α
Other Names	TNF; TNFA; TNFSF2; Tumor necrosis factor; Cachectin; TNF-alpha; Tumor necrosis factor ligand superfamily member 2; TNF-a
Accession No.	Swiss-Prot: P01375NCBI Gene ID: 7124
Uniprot	P01375
GeneID	7124;
SDS-PAGE MW	16kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/10000

Images



Western Blot analysis of Jurkat cells using TNF- α Polyclonal Antibody

Note: This product is for in vitro research use only