

## TNFAIP3 Rabbit mAb

Catalog No: #49053

Package Size: #49053-1 50ul #49053-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TNFAIP3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	SN07-31
Purification	ProA affinity purified
Applications	WB, ICC/IF, IHC, FC
Species Reactivity	Hu
Immunogen Description	recombinant protein
Other Names	A20 antibody AISBL antibody MGC104522 antibody MGC138687 antibody MGC138688 antibody OTU domain containing protein 7C antibody OTU domain-containing protein 7C antibody OTUD7C antibody Putative DNA binding protein A20 antibody Putative DNA-binding protein A20 antibody TNAP3_HUMAN antibody TNF alpha-induced protein 3 antibody TNFA1P2 antibody TNFAIP 3 antibody TNFAIP3 (A20) antibody TNFAIP3 antibody Tumor necrosis factor alpha induced protein 3 antibody Tumor necrosis factor alpha-induced protein 3 antibody Tumor necrosis factor induced protein 3 antibody Tumor necrosis factor inducible protein A20 antibody tumor necrosis factor, alpha-induced protein 3 antibody Zinc finger protein A20 antibody
Accession No.	Swiss-Prot#:P21580
Uniprot	P21580
GeneID	7128;
Calculated MW	90 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C


## Application Details

WB: 1:1,000 IHC: 1:50-1:200 ICC: 1:100-1:500 FC: 1:50-1:100


## Images




Western blot analysis of TNFAIP3 on Jurkat cells lysates using anti-TNFAIP3 antibody at 1/1,000 dilution.




Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-TNFAIP3 antibody. Counter stained with hematoxylin.




ICC staining TNFAIP3 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining TNFAIP3 in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining TNFAIP3 in HepG2 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Flow cytometric analysis of HepG2 cells with TNFAIP3 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody

## Background

---

A20 is a Cys2/Cys2 zinc finger protein that is induced by a variety of inflammatory stimuli and regulates gene expression. Specifically, A20 is induced by tumor necrosis factor (TNF) and interleukin 1 (IL-1), and acts as a negative regulator of nuclear factor  $\kappa$  B (NF $\kappa$ B) gene expression. By inhibiting NF $\kappa$ B activation, A20 plays a critical role in terminating NF $\kappa$ B responses to various stimuli. Although the C-terminal region of A20 contains seven zinc finger domains, only four of these domains are required for in vitro inhibition of TNF-induced NF $\kappa$ B activation. A20 also interacts with several other proteins, such as TRAF2, TRAF6 and I $\kappa$ B kinase (IKK)  $\gamma$  protein, and can thereby inhibit cell death. TXBP151, a novel A20-binding protein, may mediate the anti-apoptotic activity of A20. Involved in the negative feedback regulation of signal transduction, A20 and A20-binding proteins may be useful as novel therapeutic tools in the treatment of a variety of diseases.

## References

---

---

Note: This product is for in vitro research use only