

## Ubr1 Polyclonal Antibody

Catalog No: #41521

Package Size: #41521-1 50ul #41521-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ubr1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Ubr1 Polyclonal Antibody detects endogenous levels of Ubr1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ubr1.
Target Name	Ubr1
Other Names	UBR1; E3 ubiquitin-protein ligase UBR1; N-recogin-1; Ubiquitin-protein ligase E3-alpha-1; Ubiquitin-protein ligase E3-alpha-I
Accession No.	Swiss-Prot: Q8IWW7NCBI Gene ID: 197131
Uniprot	Q8IWW7
GeneID	197131;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

## Images

Western Blot analysis of Jurkat cells using Ubr1 Polyclonal Antibody



Note: This product is for in vitro research use only