

USP19 Polyclonal Antibody

Catalog No: #41524

Package Size: #41524-1 50ul #41524-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	USP19 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	USP19 Polyclonal Antibody detects endogenous levels of USP19 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human USP19.
Target Name	USP19
Other Names	USP19; KIAA0891; ZMYND9; Ubiquitin carboxyl-terminal hydrolase 19; Deubiquitinating enzyme 19; Ubiquitin thioesterase 19; Ubiquitin-specific-processing protease 19; Zinc finger MYND domain-containing protein 9
Accession No.	Swiss-Prot: O94966NCBI Gene ID: 10869
Uniprot	O94966
GeneID	10869;
SDS-PAGE MW	130kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using USP19 Polyclonal Antibody

Note: This product is for in vitro research use only