

USP44 Polyclonal Antibody

Catalog No: #41526

Package Size: #41526-1 50ul #41526-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USP44 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	USP44 Polyclonal Antibody detects endogenous levels of USP44 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human USP44.
Target Name	USP44
Other Names	USP44; Ubiquitin carboxyl-terminal hydrolase 44; Deubiquitinating enzyme 44; Ubiquitin thioesterase 44; Ubiquitin-specific-processing protease 44
Accession No.	Swiss-Prot: Q9H0E7NCBI Gene ID: 84101
Uniprot	Q9H0E7
GeneID	84101;
SDS-PAGE MW	81kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of 293T cells using USP44 Polyclonal Antibody

Note: This product is for in vitro research use only