

Vinculin Polyclonal Antibody

Catalog No: #41534

Package Size: #41534-1 50ul #41534-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vinculin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Vinculin Polyclonal Antibody detects endogenous levels of Vinculin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Vinculin around the non-phosphorylation site of Y821.
Target Name	Vinculin
Other Names	VCL; Vinculin; Metavinculin
Accession No.	Swiss-Prot: P18206NCBI Gene ID: 7414
Uniprot	P18206
GeneID	7414;
SDS-PAGE MW	123kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Vinculin Polyclonal Antibody

Note: This product is for in vitro research use only