

YAP Polyclonal Antibody

Catalog No: #41540

Package Size: #41540-1 50ul #41540-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	YAP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	YAP Polyclonal Antibody detects endogenous levels of YAP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human YAP.
Target Name	YAP
Other Names	YAP1; YAP65; Yorkie homolog; 65 kDa Yes-associated protein; YAP65
Accession No.	Swiss-Prot: P46937NCBI Gene ID: 10413
Uniprot	P46937
GeneID	10413;
SDS-PAGE MW	67kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using YAP Polyclonal Antibody



Western Blot analysis of NIH-3T3 HELA HepG2 COLO205 cells using YAP Polyclonal Antibody

Note: This product is for in vitro research use only